

Rabies Submission Form

- 1. Use one form for each specimen tested.
- 2. Complete all the information requested.
- 3. Label each specimen to correspond with the submission form.
- 4. Please read the submitting instructions given on the reverse side.

Laboratory Number			
Date Received			
Time Received			
Method of Transport			
For laboratory use only			

Animal Information				
Animal species Breed				
Description: Color				
Died ☐ Was killed ☐ Date of collection				
Did the domestic animal have a current rabies vaccination? YES □ NO □ N/A □				
Owner of suspect animal Phone				
ddress County				
() of () total specimens submitted				
NO EXPOSURE □				
Exposed Information Date of Exposure				
Name, Address and Phone Number* *(Phone number required if exposure has occurred.)	County Where Exposure Occurred	Type of Exposure	Site of Human Exposure	
Human □ Animal □		Bite Scratch Other		
Human □ Animal □		Bite Scratch Other		
Human □ Animal □		Bite ☐ Scratch ☐ Other ☐		
Human □ Animal □		Bite Scratch Other		
Submitting Agency Information				
Submitted by				
Address Phone				
Send Report to				
Address Phone				
Brief Description of Exposure/Comments (Please include where the animal was found)				
Is a fax report requested? YES \(\sigma\) NO \(\sigma\)	Is a telephone report of negative test results requested? (for emergency cases only please) YES NO			
Please give fax number here:	Please give telephone number here: After hours telephone number:			
	After hours telephone number:			

Printed by Authority of the State of Illinois P.O. #555254 500M 8/15

GUIDELINES FOR SUBMISSION OF RABIES SPECIMENS

Specimen Collection and Removal of Head and Brain

- 1. Submit only the head of ALL SKUNKS and other suspect animals.
- 2. Submit only the entire brain from large animals (e.g., horse). Bovine, deer, elk or ovine must be sent to the Department of Agriculture.
- 3. It is acceptable to submit the whole carcass from animals weighing less than 2 pounds (e.g., bat).
- 4. Live animals will not be accepted under any circumstances.
- 5 Isolate and retain the carcasses of rabies suspects until testing is completed and the report is obtained.
- 6. Head/brain removal should be performed by individuals who are immunized against rabies and trained in the proper precautions to avoid exposure, such as use of personal protection equipment.
- 7. Clean and disinfect all instruments used for head removal after each use.
- 8 For advice on testing of animals, please contact your local health department (if not available, please call 217-782-2016).
 - · Healthy domestic animals that have not exposed a person or owned animal will not be tested.
 - · Animals showing signs and symptoms of rabies should be tested, even if an exposure did not occur.
 - Submission of biting or non-biting rodents (squirrels, mice, rats, hamsters, gerbils, etc.) or rabbits is discouraged unless unusual circumstances surround the bite.
 - Animals that have bitten either a person or owned animal can be tested. Healthy dogs, cats or ferrets that bite a person or owned animal after being provoked can be confined for 10 days rather than be tested.
 - Any bat that may have exposed a person or owned animal should be tested even if a bite is not known to have occurred (e.g., found in a room with an infant or sleeping, intoxicated or otherwise incapacitated person).
 - Scratches alone have not been documented to transmit rabies; if wound(s) could have been made by scratches and/or bites, and the animal is otherwise considered a candidate for testing, or if the scratches were made by a bat, submit the animal for testing.
- 9. Rabies Submission Form: Provide all information requested on the form, if available. Be certain that a complete history accompanies each specimen submitted (see reverse side). Indicate on the submission form if more than one animal is submitted. Mark individual specimens with the same identification (e.g., #3 of 5 total animals submitted).

Packaging and Transportation

- 1. Place the specimen in a primary container such as a sealed plastic bag. Bats and other small animals must be placed inside a clear zip-lock bag or equivalent. Place the primary container in a watertight secondary container with sufficient absorbent material to assure complete absorption if leakage should occur. If plastic bags are used, please TRIPLE bag each specimen. Place the test requisition on the outside of the secondary container. Package forms so that they remain clean.
- 2. Place the secondary container inside a rigid outer container with sufficient cold packs to maintain refrigeration temperatures. The minimum volume of cold packs surrounding the specimen is double the size of the specimen (never use ice cubes).
- 3. Ample insulation surrounding the specimen and cold packs is essential (e.g., styrofoam packing material or newspaper).
- 4. **Do not freeze specimens** unless the specimen is held over a long weekend (weekend + holiday) or if the animal is decomposed. Cold, not frozen is the general rule. Do not package rabies specimens with dry ice since dry ice may freeze the specimens.
- 5. Pack rabies specimens separately from other types of laboratory specimens (even if they are going to the same lab).
- 6. **Messenger/Courier by ground transport** Place the secondary container and test requisition in the shipping container. The shipping container must be rigid such as a cooler and labeled with the **UN 3373** Biological Substance Category B marking. Close securely. **Styrofoam coolers must be placed in a cardboard box.**
- 7. Commercial carrier by ground/air transport Place the secondary container and test requisition in the shipping container.

 Label the outer shipping container with the appropriate Illinois Department of Public Health laboratory address. Complete the return address section to include the name of the person shipping the package, business name and address and a business phone number. The shipping container must include the
 - UN3373 Biological Substance Category B marking.
- Do not ship specimens on Fridays or the day before a state holiday unless prior approval is obtained through

You as the shipper, not the transport company, are responsible for the shipment until the package reaches its destination.

the IDPH Communicable Disease control section. Do not ship on Thursdays unless sent by guaranteed, overnight delivery. Specimens held by a commercial carrier over a weekend are frequently received in various stages of decomposition that may preclude the possibility of reliable test results. In case of emergency, please notify the appropriate laboratory during working hours, 8:30 a.m. to 4 p.m, Monday through Friday. For an after-hours emergency, notify IEMA at 800-782-7860. Ask to speak with the Illinois Department of Public Health Duty Officer.

- 9. Send specimens by the most expeditious means available, preferably by direct messenger. However, if a commercial carrier is used, the shipper is responsible for conforming to any special packaging requirements of the carrier.
- 10. The IDPH has no pickup and delivery service and assumes no responsibility for any specimen until it is received and identified.

ILLINOIS DEPARTMENT OF PUBLIC HEALTH LABORATORIES

Illinois Department of Public Health Division of Laboratories 1155 S. Oakland Ave. P.O. Box 2797 Carbondale, IL 62901 618-457-5131 Illinois Department of Public Health Division of Laboratories 2121 W. Taylor St. Chicago, IL 60612 312-793-4760 Illinois Department of Public Health Division of Laboratories 825 N. Rutledge St. P.O. Box 19435 Springfield, IL 62702 217-782-6562

Animal testing other than for rabies is **not** performed by the Illinois Department of Public Health laboratories.

ILLINOIS DEPARTMENT OF AGRICULTURE LABORATORIES

Animal Disease Laboratory 2100 S. Lake Storey Road Galesburg, IL 61401 309-344-2451